3M ESPE	
<div> <div></div> <div>Ketac™ Bond Apicap™</div> </div>	
 Glass Ionomer Base Material	
 Glasionomer Unterfüllungsmaterial	
 Ciment verre-ionomère pour fond de cavité	
 Materiale da sottofondo vetro ionomero	
 Material de obturación de ionómero de vidrio	
 Material ionómero vítreo para obturação de base	
 Glasionomeer onderlaagcement	
 Υαλιονομερές υλικό για ουδέτερα στρώματα	
 Glasjonomert underfyllningsmaterial	
 Lasi-ionomeeri-alustäytämateriaali	
 Glasionomer-underfyldningsmateriale	
 Glassionomer underföringsmaterial	
<p>Instructions for Use</p> <p>Gebrauchsinformation</p> <p>Mode d'emploi</p> <p>Informazioni per l'uso</p> <p>Información de uso</p> <p>Instruções de uso</p>	<p>Gebruiksinformatie</p> <p>Πληροφορίες χρήσεως</p> <p>Bruksanvisning</p> <p>Käyttöinfo</p> <p>Brugsanvisning</p> <p>Bruksinformasjon</p>
<div> <div></div> <div>3M ESPE AG</div> <div>Dental Products</div> <div>D-82229 Seefeld – Germany</div> </div>	
<div> <div></div> <div>3M ESPE</div> <div>Dental Products</div> <div>St. Paul, MN 55144-1000 U.S.A.</div> </div>	
3M ESPE Customer Care/MSDS Information: U.S.A. 1-800-634-2249 and Canada 1-888-363-3685.	
3M, ESPE, Apicap, CapMix, Ketac and RotoMix are trademarks of 3M or 3M ESPE AG.	
© 3M 2010. All rights reserved.	

ENGLISH

Product Description

Ketac Bond Apicap, manufactured by 3M ESPE, is a fast-setting glass ionomer cement supplied in capsules and has a dense consistency. The feature of the material to chemically bond to enamel and dentin provides for gentle preparation. Ketac Bond Apicap releases fluoride ions and is radio-opaque. Ketac Bond Apicap is available in one color, yellow. Each capsule provides 0.1 ml of usable material.

For details on all products mentioned below please refer to the corresponding Instructions for Use.

Instructions for Use should not be discarded for the duration of the product use.

Indications

- Base linings under composite and amalgam
- Extended fissure sealing without acid etching
- Core build-ups

Precautionary Measures

3M ESPE MSDSs can be obtained at www.mmm.com or from your local subsidiary.

Preparation

- Only remove the carious substance; no undercuts required.

Protection of the Pulp

- Cover areas in close proximity to the pulp by applying spots of a hard-setting calcium hydroxide material, in order to protect the pulp from irritation.

Conditioning

- For optimal chemical bonding to the teeth the smear layer produced in the preparatory steps must be carefully removed by applying Ketac Conditioner, manufactured by 3M ESPE, to the prepared surfaces and allowing the substance to react for 10 sec.
- Then rinse with copious amounts of water.
- Blow the cavity dry in only 2–3 short intervals with water and oil-free air or dab dry with cotton pellets. Do not over-dry! **The cavity should be just so dry that the surface has a matt shiny appearance.** Excessive drying can result in post-operative sensitivity after filling.
- Avoid further contamination.

Activation

- Place the Apicap™ Aktivator, manufactured for 3M ESPE, on a sturdy surface and insert the Apicap capsule in the activator.
- Using the ball of your hands, depress the activator lever firmly and **completely to the stop** and hold it down for 2 to 4 seconds.
 - Firmly and completely depressing the lever to the stop and holding it in this position is the only way to ensure that all of the liquid is pressed into the powder. If inadequate pressure is applied or if the lever is not depressed completely it is possible that too little liquid will be mixed with the powder, increasing the viscosity and changing the product characteristics.

Mixing

- Mix the capsules at approx. 4,300 rpm, in a high frequency mixing device such as CapMix™ (for 10 sec) or the RotoMix™ rotating mixer (for 8 sec), and apply with the Apicap™ Applier, all products manufactured by or for 3M ESPE. Both the Apicap Aktivator and the Apicap Applier are color-coded in red.

Caution

- Delays between the steps of activation, mixing, and application must be carefully avoided, since the application of the material may be impeded or prevented by the material beginning to set.

Application

- Process Ketac Bond Apicap at room temperature (20–25° C/68–77° F).
- Ensure that there is no contamination from water, saliva, etc., entering the working area for the entire duration of application; the working area must be kept dry.
- Aesthetic fillings can be produced with the so-called sandwich technique, in which the glass ionomer cement base filling replacing the dentin is covered with a layer of composite material according to the following procedure:
- Apply the desired quantity of Ketac Bond Apicap with a suitable instrument leaving the margins of the dentin exposed. After the Ketac Bond Apicap material has set, etch both enamel and lining and then apply the desired composite using a suitable bonding material. Process this material in accordance with manufacturer recommendations. Lining located below amalgam does not have to be etched and banded.
- For fissure sealing and core build-ups, fill the cavity using a suitable instrumt. Wait a minimum of 4 minutes from start of mixing prior to commencing preparation of the hardened Ketac Bond Apicap with rotary instruments.
- Ketac Bond Apicap adheres to metal instruments. Therefore, residual quantities of the cement should be washed off with cold water before they set on the instrument.

Times

At an ambient temperature of 23° C/74° F and 50% relative humidity the following times apply:

Mixing (RotoMix)	min:sec
	00:08
or high frequency mixer, e.g. CapMix	00:10
Processing incl. mixing	01:45
Setting from start of mixing	04:00
Setting inside the mouth	02:15

At higher temperatures than indicated above the time available for processing is shortened, whereas it is prolonged at lower temperatures.

Storage and Shelf-life

Store the product at 15–25° C/59–77° F. Do not use after the expiration date.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manu- facture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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